



Socky and Butler dressing aids

Socky stocking aids make it possible for persons with impaired mobility to put on socks, stockings and support stockings. This range is particularly suitable for users with back, hip or knee problems.

The Butler buttoners help users with impaired mobility in the arms or hands to manage buttons.



Socky stocking aid, short

Socky short stocking aid, can be used to put on long or short stockings or socks. It is small and flexible and therefore easy to bring. The stocking slides over the nylon fabric, which is then easily pulled over the foot.



Socky stocking aid, short

Description	Item no.	Material:	Fabric: Polyamide
Socky stocking aid, short	8060 1002	Plastic parts: Polypropylene	
		Weight:	94 g
		Length:	40 cm
		Design:	HM Center
		Patent:	Socky short is patented

Socky stocking aid, long

Socky long stocking aid, makes it easy for users to put on long or short stockings, socks or tights. For tights, use two Socky stocking aids connected with the linking device. The plastic sheet slides easily on all floors and has a grip which makes it easier to pull up the stocking.

Description	Item no.	Material:	Fabric: Polyamide
Socky stocking aid, long	8060 1001	Plastic parts: Polypropylene	
Linking device	8060 1021	Weight:	146 g
		Length:	60 cm
		Design:	HM Center
		Patent:	Socky long is patented



Socky stocking aid, long



Socky stocking aid for support stocking and compression stockings

This stocking aid is designed for putting on support stockings of various lengths and types. It has large handles which facilitates also for an assisting carer.



Socky stocking aid, for support and compression stocking

Description	Item no.	Material:	Polyamide fabric
Socky stocking aid for support and compression stockings	8060 1003	Weight:	90 g
		Length:	100 cm
		Design:	HM Center

Butler, available in three variants

The Butler buttoners are for people who cannot manage buttons due to impaired mobility in the arms or hands.

Three types of buttoners are available: a buttoner with spherical grip for small buttonholes, and two buttoners with combination loops for all other buttonholes. One has a long handle. The other has a spherical grip.



Butler combi buttoner



Butler buttoner



Butler buttoner combi long



Dimensions:	Butler small and Butler combi: 13 cm
	Butler combi long: 22 cm
Material:	Handle: Polystyrene
	Shank: Epoxy-coated metal or galvanised wire
Weight:	Butler small and Butler combi: 25 g
	Butler combi long: 18 g
Design:	HM Center

Description	Item no.
Butler buttoner for small buttonholes	8060 2001
Butler combi buttoner for all other buttonholes	8060 2003
Butler buttoner combi long for all other buttonholes	8060 2004